

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	02-352-6000	Revision: 4
Remarks		single use only

Specification

Product	Test tube PP with fitted screw cap PE, 12ml, Ø16x100 mm, round bottom, tube transparent, cap natural, grad., sterile R, max. RCF 4.000g, Autocl. up to 121°C (no cap)
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height (without screw cap): 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	grad.
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	tube: -196°C up to +121°C
	cap: -50°C up to +80°C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	