

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	02-421-0000	Revision: 2
Remarks	Fitting plug: 02-673-0000	single use only

Specification

Product	Test tube PS, conical bottom, 12ml, Ø16x105 mm, transparent, max. RCF 3.000g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height: 105,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	12ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000
	pcs per carton: 2.000
	pcs per pallet: 40.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	