

Technical Datasheet

Date: 06.11.24

Page 1 of 1

EN ISO 13485 certified		
Article - No:	02-552-6000	Revision: 3
Remarks		single use only

Specification

Product	Test tube PP with fitted screw cap PE, 30ml, Ø25x107 mm, conical bottom, skirted, tube transparent, cap natural, grad., sterile R, max. RCF 5.000g, autocl. up to 121°C
Dimensions	outer diameter: 25,00 mm (+/- 0,2 mm)
	height: 107,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	grad.
Nominal Capacity	30ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	tube: -196°C up to +121°C
	cap: -50°C up to +80°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 5.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months
Packaging	pcs per bag: 100
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	