

Technical Datasheet

Date: 13.10.23

Page 1 of 2

EN ISO 13485 certified		
Article - No:	02-572-4001	Revision: 8
Remarks		single use only

Specification

Product	Centrifuge tube PP with screw cap PE, 50ml, Ø30x115 mm, conical, skirted, tube transp., cap green, grad., writ. area, np pcr ready, max. RCF 7.000g, Autocl. up to 121°C (no cap)
Dimensions	height: 115,00 mm (+/- 0,2 mm)
	diameter: 30,00 mm (+/- 0,2 mm)
Colour	tube transp., cap green
Graduation	grad., writ. area
Nominal Capacity	50ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	cold resistant - 80 °C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 7.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	free of detectable phthalate
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 25
	pcs per carton: 250
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C

Technical Datasheet

Date: 13.10.23

Page 2 of 2

	pal. not stackable
Approved by nerbe plus quality department	