

Technical Datasheet

Date: 16.08.24

Page 1 of 2

EN ISO 13485 certified		
Article - No:	02-572-9001	Revision: 10
Remarks		single use only

Specification

Product	Centrifuge tube PP with screw cap PE, 50ml, Ø30x115 mm, conical, skirted, cap green, grad., green, grad., writing area, np pcr ready, steril R max RCF 7.000g, Autocl. up to 121°C (no cap)
Dimensions	height: 115,00 mm (+/- 0,2 mm)
	diameter: 30,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap green
Graduation	grad., writ. area
Nominal Capacity	50ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	cold resistant - 80 °C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 7.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	free of detectable phthalate
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 25
	pcs per carton: 250

Technical Datasheet

Date: 16.08.24

Page 2 of 2

	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	