

Technical Datasheet

Date: 10.12.24

Page 1 of 2

EN ISO 13485 certified		
Article - No:	04-083-0096	Revision: 6
Remarks	for Roche Lightcycler 480	single use only

Specification

Product	PCR-plate PP, 96x0,2ml, half-skirted, low profile white, np pcr ready, max. RCF 3.000g, autocl. up to 121°C
Dimensions	ANSI/SLAS 1 & 4
	length: 127,76 mm (+/- 0,5 mm)
	width: 85,48 mm (+/- 0,5 mm)
Colour	white
Nominal Capacity	96x0,2ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Packaging	pcs per inner carton: 10
	pcs per carton: 100
	pcs per pallet: 4.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 10.12.24

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

