

Technical Datasheet

Date: 06.11.24

Page 1 of 2

EN ISO 13485 certified		
Article - No:	04-252-9000	Revision: 4
Remarks	piercable cap	single use only

Specification

Product	Microcentrifuge tube PP, 5ml, attached cap, transparent,, grad., writ. area, np pcr ready, sterile R, max. RCF 25.000g, autocl. up to 121°C
Dimensions	outer diameter: 14,70 mm (+/- 0,2 mm)
	inner diameter: 12,60 mm (+/- 0,2 mm)
	height: 58,20 mm (+/- 0,2 mm)
Colour	transparent,
Graduation	grad., writ. area
Nominal Capacity	5ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 25.000g
Product Characteristics	biologically inert
	physiologically harmless
	constantly usable at temperature of 80 °C up to 100 °C
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months
Packaging	pcs per bag: 200
	pcs per carton: 4.000
	pcs per pallet: 56.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C

Technical Datasheet

Date: 06.11.24

Page 2 of 2

	pal. not stackable
Approved by nerbe plus quality department	