

# Technical Datasheet

Date: 11.06.24

Page 1 of 1

EN ISO 13485 certified		
Article - No:	05-500-0015	Revision: 4
Remarks		single use only

## Specification

Product	Micro screw cap tube PP, 1,5ml, with enclosed screw cap PP with O-ring, tube natural, cap transparent, max. RCF 20.000g
Colour	tube natural, cap transparent
Nominal Capacity	1,5ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Centrifugation	max. RCF 20.000g
Purity Grade	free of RNase and DNase free of pyrogens
Packaging	pcs per bag: 500 pcs per carton: 4.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	