

# Technical Datasheet

Date: 27.10.23

Page 1 of 1

|                        |             |                 |
|------------------------|-------------|-----------------|
| EN ISO 13485 certified |             |                 |
| Article - No:          | 05-900-1010 | Revision: 3     |
| Remarks                |             | single use only |

## Specification

|                                           |                                                                                                                             |
|-------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product                                   | Cryo vial PP, screw cap PE, 1,0ml, self-standing, transparent, grad., sterile R, autocl.                                    |
| Dimensions                                | height: 36,1mm                                                                                                              |
| Colour                                    | transparent                                                                                                                 |
| Graduation                                | grad.                                                                                                                       |
| Nominal Capacity                          | 1,0ml                                                                                                                       |
| Material Characteristics                  | tube: polypropylene (PP)<br>cap: polyethylene (PE)                                                                          |
| Temperature range                         | -196°C - 121°C                                                                                                              |
| Autoclavability                           | autocl.                                                                                                                     |
| Purity Grade                              | free of PCR inhibitors<br>free of human DNA<br>free of RNase and DNase<br>free of pyrogens                                  |
| Sterilisation Process                     | radiation sterilization SAL 10-6                                                                                            |
| Packaging                                 | pcs per bag: 100<br>pcs per carton: 1.000<br>pcs per pallet: 99.000                                                         |
| Storage                                   | dry condition: 40%-60% rel. humidity<br>bags / cases orig. closed<br>storage temp. +2°C to max. +30°C<br>pal. not stackable |
| Approved by nerbe plus quality department |                                                                                                                             |