

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	05-900-1050	Revision: 3
Remarks		single use only

Specification

Product	Cryo vial PP, screw cap PE, 5ml, self-standing, transparent, grad., sterile R, autocl.
Dimensions	height: 88,5mm
Colour	transparent
Graduation	grad.
Nominal Capacity	5ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-196°C - 121°C
Autoclavability	autocl.
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 100
	pcs per carton: 500
	pcs per pallet: 49.500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	