

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	05-902-0009	Revision: 5
Remarks		single use only

Specification

Product	Insert for cryo tube lid, 5 different colors, autocl.
Colour	5 different colors
Temperature range	-196°C - 121°C
Autoclavability	autocl.
Packaging	pcs per bag: 500
	pcs per carton: 500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	