

Technical Datasheet

Date: 13.10.23

Page 1 of 2

EN ISO 13485 certified		
Article - No:	06-372-2018	Revision: 12
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, refill system, 0,1-20µl, super slim, 96 pcs/tray, transparent tips, translucent gray tray, grad., np premium bio, sterile R, autocl. up to 121°C
Dimensions	length: 45,8 mm (+/- 0,2 mm)
	inner diameter: 0,38 mm (+/- 0,2 mm)
Colour	transparent tips, translucent gray tray
Graduation	grad.
Nominal Capacity	0,1-20µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens and endotoxins
	Sterile R
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960
	pcs per carton: 9.600
	pcs per pallet: 518.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 13.10.23

Page 2 of 2

Approved by nerbe plus quality department