

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-895-0200	Revision: 5
Remarks	suitable for: Eppendorf and Biohit	single use only

Specification

Product	Pipette tip PP "classic", 5.000µl, universal fit (Eppendorf, Biohit, Socorex), transparent, grad. np pcr ready, autocl. up to 121°C
Dimensions	length: 120,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	5.000µl
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 250
	pcs per carton: 2.500
	pcs per pallet: 50.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

