

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	07-962-0301	Revision: 7
Remarks		single use only

Specification

Product	Gel-loading pipette tip PP, universal fit, 200µl, round tip, 96 pcs/rack, np pcr ready
Dimensions	length: 67,20 mm (+/- 0,2 mm)
Nominal Capacity	200µl
Material Characteristics	polypropylene (PP)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per inner cs: 768
	pcs per carton: 3.072
	pcs per pallet: 122.880
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	