

Technical Datasheet

Date: 22.07.24

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-452-6000	Revision: 11
Remarks	max. volume: 140ml / working volume: 125 ml	single use only

Specification

Product	Screw cap container PP, 125ml, Ø58X80mm, cap pre-assembled, container transparent, cap green, grading 3 sided, every 10ml, sterile EO
Colour	container transparent, cap green
Graduation	grading 3 sided, every 10ml
Nominal Capacity	125ml
Material Characteristics	container: polypropylene (PP) cap: high-density polyethylene (HDPE)
Temperature range	-40°C - 120°C
Product Characteristics	leak proof of not less than 95 KPa (0,95 bar, 14 psi) biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	RNase, DNase, Pyrogene free
Sterilisation Process	gas sterilisation by ethyleneoxid, SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 150 pcs per carton: 150 pcs per pallet: 4.800
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	