

Technical Datasheet

Date: 13.10.23

Page 1 of 2

EN ISO 13485 certified		
Article - No:	11-702-1010	Revision: 1
Remarks	suitable for formalin and chloroform	single use only

Specification

Product	Screw cap container PP, 500ml, Ø90x91 mm, with fitted screw cap PP, with PE-pad, white, autoclavable up to 121°C (without pad)
Dimensions	container:
	diameter: 89,20 mm (+/- 0,2 mm)
	height: 90,50 mm (+/- 0,2 mm)
	cap:
	diameter: 96,00 mm (+/- 0,2 mm)
	height: 15,70 mm (+/- 0,2 mm)
	total height container and cap: 91,00 mm (+/- 0,2 mm)
Colour	white
Nominal Capacity	500ml
Material Characteristics	container: polypropylene (PP)
	cap: polypropylene (PP)
	pad: polyethylene (PE)-foam (400 kg/cbm)
Temperature range	-80°C - 70°C
Autoclavability	autocl. up to 121°C (without pad)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per cs: 72
	pcs per carton: 72
	pcs per pallet: 1.728
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed

Technical Datasheet

Date: 13.10.23

Page 2 of 2

	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	