

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-712-1010	Revision: 2
Remarks	suitable for formalin and chloroform	single use only

Specification

Product	Screw cap container PP, 1.250ml, Ø106x159 mm, with fitted screw cap PE, with EPE-pad, white, autoclavable up to 121°C (without pad)
Colour	white
Nominal Capacity	1.250ml
Material Characteristics	container: polypropylene (PP)
	cap: polypropylene (PP)
	pad: expanded polyethylene (EPE)
Temperature range	-80°C - 70°C
Autoclavability	autocl. up to 121°C (without pad)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
Packaging	pcs per cs: 35
	pcs per carton: 35
	pcs per pallet: 700
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	