

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-842-1100	Revision: 1
Remarks		single use only

Specification

Product	Multi purpose container PP, 1.000ml, Ø110x132 mm, with enclosed snap-on lid PE , transparent, graduated
Dimensions	outer diameter: 110,00 mm (+/- 0,2 mm)
	inner diameter: 107,00 mm (+/- 0,2 mm)
	height (with lid): 132,00 mm (+/- 0,2 mm)
	height (without lid): 131,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	1.000ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-80°C - 75°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Packaging	pcs per bag: 250
	pcs per carton: 250
	pcs per pallet: 3.750
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	