

Technical Datasheet

Date: 13.10.23

Page 1 of 2

EN ISO 13485 certified		
Article - No:	11-902-2100	Revision: 1
Remarks		single use only

Specification

Product	Hinged lid container PP, 50ml, Ø33x80 mm, natural Moulded rings every 10 ml, writing area, np pcr ready, autoclavable up to 121°C
Dimensions	outer diameter (without lid): 33,00 mm (+/- 0,2 mm)
	inner diameter: 31,00 mm (+/- 0,2 mm)
	height (without lid): 78,00 mm (+/- 0,2 mm)
	height (with lid): 80,00 mm (+/- 0,2 mm)
Colour	natural
Graduation	Moulded rings every 10 ml, writing area
Nominal Capacity	50ml
Material Characteristics	polypropylene (PP)
Temperature range	-20°C - 120°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 650
	pcs per carton: 650
	pcs per pallet: 13.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 13.10.23

Page 2 of 2

Approved by nerbe plus quality department