

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-250-5000	Revision: 4
Remarks		single use only

Specification

Product	Straw pipette PP, 1ml, Ø4,5x190 mm, 25 pcs/bag, transparent, sterile R
Dimensions	diameter: 4,50 mm (+/- 0,2 mm)
	length: 190,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	1ml
Material Characteristics	polypropylene (PP)
Temperature range	-20°C - 120°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per inner cs: 1.000
	pcs per carton: 12.000
	pcs per pallet: 240.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	