

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-541-5103	Revision: 3
Remarks		single use only

Specification

Product	Spreader PS, Ø3,5x145 mm, width 35 mm, T-shaped, 5 pcs/peelpack, blue, sterile EO
Dimensions	diameter: 3,50 mm (+/- 0,2mm)
	length: 145,00 mm (+/- 0,2mm)
	width spatula: 38,00 mm (+/- 0,2mm)
Colour	blue
Material Characteristics	polystyrene (PS)
Temperature range	-10°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	gas sterilisation by ethyleneoxid, 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per inner cs: 500
	pcs per carton: 2.000
	pcs per pallet: 56.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	